DESIGN QUALIFICATION

FOR

AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE

DOCUMENT No.: FTIL- RASTAGENE -20-DQ-VD-002



Life Engineering[™]

717, Janki Centre, Off Veera Desai Road, Andheri (W),

Mumbai 400053, India Phone:-+912261592900

Website:- www.fabtechnologies.com Email:- sales@fabtechnologies.com

TABLE OF CONTENTS



M/s. RASTAGENE BIOPHARMACEUTICALS CO

| Life Engineering | | | |
|------------------|---|---------------|---|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 2 of 16 | |

| 1 | APPROVAL SHEET | | | | | | |
|-----|---------------------|-----------------------------------|---|--|--|--|--|
| 2 | SYSTEM INFORMATION4 | | | | | | |
| 1.0 | Овјест | IVE4 | ŀ | | | | |
| 2.0 | SCOPE | 4 | ŀ | | | | |
| 3.0 | REFERE | NCES | ; | | | | |
| 4.0 | RESPON | ISIBILITIES5 | ; | | | | |
| | 4.1 | CLIENT'S | 5 | | | | |
| | 4.2 | Manufacturer's5 | | | | | |
| 5.0 | System | DESCRIPTION6 | , | | | | |
| 6.0 | HARD | WARE SPECIFICATION7 | 7 | | | | |
| | 6.1 | MECHANICAL HARDWARE SPECIFICATION | | | | | |
| | 6.2 | ELECTRICAL SPECIFICATION | | | | | |
| 7.0 | TECHNI | CAL SPECIFICATIONS10 |) | | | | |
| 8 | SAFETY | FEATURES | } | | | | |
| 9 | | Y SPECIFICATIONS14 | | | | | |
| 10 | REPORT | T APPROVAL15 | ; | | | | |
| 11 | ABBREV | /IATION16 | Ś | | | | |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS (| CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 3 of 16 | |

1 APPROVAL SHEET

| M/s. FABTECH TECHNOLOGIES INTERNATIONAL LTD. | | | | | |
|--|------|------------|-------------|-----------|------|
| | Name | Department | Designation | Signature | Date |
| Prepared by | | | | | |
| Checked by | | | | | |
| Approved by | | | | | |

| | M/s. RASTAGENE BIOPHARMACEUTICALS CO. | | | | | |
|-------------|---------------------------------------|------------|-------------|-----------|------|--|
| | Name | Department | Designation | Signature | Date | |
| Reviewed by | | | | | | |
| Reviewed by | | | | | | |
| Approved by | | | | | | |

REVISION HISTORY

| SR. NO. | REVISION | DATE | REVISION SUMMARY |
|---------|----------|------|------------------|
| | | | |
| | | | |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 4 of 16 | |

2 SYSTEM INFORMATION

| EQUIPMENT | AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE |
|------------------|--|
| MANUFACTURER | FABTECH TECHNOLOGIES INTERNATIONAL LTD. |
| CUSTOMER | M/s. RASTAGENE BIOPHARMACEUTICALS CO |
| SITE | IRAN |

1.0 OBJECTIVE

This document summarizes the Design Qualification of **Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit** that will be installed at the Customer. The purpose of the Design Qualification (DQ) is to provide the specifications and operation requirements of the equipment based on which it is designed. It is used to specify the requirements for the Installation Qualification of the equipment.

2.0 SCOPE

This DQ is applicable to **Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit** which will be used for fill liquid in vials and insert rubber plugging.

The system boundary is defined as including the **Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit,** it will define the responsibilities, acceptance criteria, basis of design, technical specifications and list of major bought out parts, utility requirements, safety & cGMP features and documentation requirements.

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS (| CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 5 of 16 | |

3.0 REFERENCES

User requirement specification (URS) given by customer.

4.0 RESPONSIBILITIES

4.1 CLIENT'S

- ➤ For providing complete technical specifications for design along with Purchase Order.
- > To qualify equipment along with manufacturer.

4.2 MANUFACTURER'S

- ➤ To manufacture and supply the equipment incorporating all the specifications as per purchase order and user requirement specification (URS).
- ➤ To assist client for successful installation & commissioning at the site.
- ➤ To design, engineer and provide the complete technical details of the equipment pertaining to its design qualification viz.
- Machine overview
- > Equipment orientation with layout
- > Specification of the sub-components/ bought out items, their make, model, quantity.
- And backup records/brochures.
- Details of Utilities.
- ➤ Identification of components for calibration.
- Material of construction of all components.
- Brief process description.
- Safety features and alarms.
- Pre-installation requirements.

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 6 of 16 | |

5.0 System Description

The Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit is in the filling Room of the Injectable Facility with restricted access. The Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit is located such that it can be attended easily for routine operations, monitoring, and sampling and maintenance functions.

The Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit is designed to fill sterile liquid Injectable with different volume in different sizes of Vial and rubber stoppered the same as well pressing of the Rubber Stoppered Vial.

The speed of the machine is 40 vials per minute depending upon the size/shape of the Vials.

The Automatic Isolator Compatible Two Head Injectable Liquid Filling Machine with Rotary Disc Type Rubber Stoppering Unit works on volumetric Filling principle giving guarantee of high volumetric accuracy with minimal spillage.

Two inlet nozzles connected to the filling vessel via silicone tubing. Outlet of Peristaltic Pump connected to filling nozzle via silicon tube. Peristaltic Pump are suck the liquid from vessel and pump out into empty vial via nozzle. Nozzle will include pre and post nitrogen facility.

As soon as empty vial comes near to filling station, motion star wheel will stop the vial and start filling. As soon as filling cycle completes, motion star wheel will release the filled vial and transport conveyor for rubber stoppering process.

Now filled vial pick-up by intermittent motion star wheel and bring the filled vial at rubber stoppering station. Rubber stoppering unit is rotary disc type (half and full bunging). Roller will firmly stoppering rubber on vial in motion. Once this cycle complete, vial comes on the conveyor again to release. And collected in out feed turn table.

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 7 of 16 | |

6.0 HARDWARE SPECIFICATION

Note: The identifiers/ Tag numbers for the major components will be listed during the qualification of the equipment.


6.1 MECHANICAL HARDWARE SPECIFICATION

| SR. NO. | ITEM NAME | мос |
|------------|----------------------------|--|
| 1. | MACHINE FRAME STRUCTURE | S.S.304 PIPE (50 X 50) |
| 2. | MACHINE TOP PLATE | SS 316 L (20 MM THICKNESS) |
| 3. | DOOR SET | S.S.304 |
| 4. | STAR WHEEL | DELRIN |
| 5. | INFEED GUIDE | DELRIN |
| 6. | NOZZLE STAND | DELRIN |
| 7. | ALL COVERS | S.S.316L WITH SILICON SHEET FOR WASHING COMPATIBLE |
| 8. | TURN TABLE PLATE | DELRIN |
| 9. | CHANGE PARTS | DELRIN |
| 10. | SYRINGES | SS 316 L |
| 11. | BOWL | SS 316 L |
| 12. | CHUTE | SS 316 L |
| 13. | NEEDLES | SS 316 L |
| 14. | CONVEYOR | DELRIN |
| 15. | DRAINAGE PIPES | SS 304 |
| 16. | BUFFER TANK 3 LITER | SS 316L |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | CO |
|--|---|---------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 8 of 16 | |

6.2 ELECTRICAL SPECIFICATION

| SR. NO. | ITEM NAME | MAKE | MODEL NO |
|------------|----------------------------|------------------|---|
| 1. | IN FEED TURN TABLE GEARBOX | BONFIGLIOLI | W-63 UFC1 100 P71 B5B3 |
| 2. | IN FEED TURN TABLE MOTOR | BONFIGLIOLI | BE 71B 4 230/400-50 IPSS CLF B5 BIS |
| 3. | OUTFEED TURN TABLE MOTOR | BONFIGLIOLI | W-63 UFC1 100 P71 B5B3 |
| 4. | OUTFEED TURN TABLE GEARBOX | BONFIGLIOLI | BE 71B 4 230/400-50 IPSS CLF B5 BIS |
| 5. | CONVEYOR MOTOR | ROTOMAQ | 120W R15 FHP GEAR MOTOR |
| 6. | SERVO | SIEMENS | V90 SERIES |
| 7. | PERISTALTIC PUMP | WATSON MARLOW | 313D PUMPHEAD, 3 ROLLERS 1.6MM WT TUBING |
| 8. | AC DRIVE | SIEMENS | 6SL3210-5BB13-7UV V20 0.37 KW |
| 9. | PLC | SIEMENS | S7 1200 |
| 10. | HMI | SIEMENS | TP - 900 COMFORT |
| 11. | MCB | SIEMENS | |
| 12. | PNEUMATIC PISTON | SMC | CQ2A20-10DM, CQ2A50-50DZ |
| 13. | PNEUMATIC 5/2 VALVE | SMC | SY7120-SL-02 |
| 14. | NITROGEN FLUSHING VALVE | SMC | VX20AAXB |
| 15. | VACUUM PUMP | BECKER | VT 4.10-G008220(G008220) |
| 16. | DOOR PROXY | SICK | |



| 17. | PHOTOELECTRIC SENSOR | SICK | 1052444 |
|-----|-------------------------|--------|--------------|
| 18. | FIBER OPTIC SENSOR | SICK | |
| 19. | AIR PRESSURE | KELLER | 0 TO 10 BAR |
| 17. | TRASNMITTER | KEEEEK | 0 TO TO BITE |
| 20. | VACUUM PRESSURE | KELLER | -1 TO 0 BAR |
| 20. | TRASNMITTER | KLLLLK | 1 10 0 DAK |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS (| СО |
|--|---|----------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 10 of 16 | |

7.0 TECHNICAL SPECIFICATIONS

7.1 Basic Machine

| PARAMETER | Description |
|------------|---|
| BASIC UNIT | MACHINE IS SUITABLE FOR 2ML AND 100ML VIAL WITH DIAMETER 24" SCRAMBLER AND UN SCRAMBLER |

7.2 CONTAINER SIZE

| PARAMETER | DESCRIPTION |
|-------------|---|
| ROUND SHAPE | 25 mm to 48 mm (Diameter) 45 mm to 100 mm (Height) |

7.3 CAPACITY / OUTPUT

| PARAMETER | DESCRIPTION |
|---------------|---------------------------------|
| MACHINE SPEED | 35 to 40 Vials / Min (For 10ML) |

7.4 MATERIAL OF CONSTRUCTION

| PARAMETER | Description |
|-------------------|-------------|
| CONTACT PARTS | SS 316 L |
| Non-Contact Parts | SS 316L |

7.5 **OPERATING SYSTEM**

| PARAMETER | DESCRIPTION |
|------------------|-------------|
| Operating System | нмі |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS (| СО |
|--|---|----------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 11 of 16 | |

7.6 OTHER FEATURES

| PARAMETER | Description |
|---------------|--|
| Gear Box | Make Bonfiglioli gear box with motor for machine |
| User Friendly | All parts for cleaning are easily removable. |
| Machine Frame | S.S.304 Pipe Structure |

7.7 SLAT CONVEYOR

| PARAMETER | Description |
|-----------|-------------------------|
| Conveyor | TO TRANSPORT THE OBJECT |

7.8 FILLER STAR WHEEL

| PARAMETER | Description |
|-------------------|--|
| FILLER STAR WHEEL | TO STOP THE VIAL ON CONVEYOR FOR FILLING PROCESS |

7.9 'FILLING HEAD

| PARAMETER DESCRIPTION | |
|-----------------------|-----------------------------------|
| FILLING HEAD | Drive mechanism to drive syringes |

7.10 MAIN DRIVE

| PARAMETER | Description |
|------------|--------------------------------------|
| Main Drive | To drive the syringes and star wheel |

7.11 Rubber Stopper Hopper & Chute

| PARAMETER | Description |
|----------------------------------|--|
| Rubber Stopper Hopper & Chute | To load bulk rubber stopper & transfer chute |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS (| CO |
|--|---|----------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 12 of 16 | |

7.12 SENSOR

| PARAMETER | Description | |
|------------|----------------------------------|--|
| Sensor – 1 | No Vial – No Fill | |
| Proxy - 1 | Door Interlock | |
| Sensor – 2 | No rubber Stopper – Machine Stop | |
| Sensor – 3 | Out Feed Vial Jam – Machine Stop | |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS CO | |
|--|---|----------------|---|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 13 of 16 | |

8 SAFETY FEATURES

The equipment shall be provided with safety features as listed below.

| SR.NO | FEATURE | SPECIFIED |
|-------|---------------------|--------------|
| 1. | Emergency Stop | Machine Stop |
| 2. | Motor Overload | Machine Stop |
| 3. | Door Open | Machine Stop |
| 4. | No vial at infeed | Machine Stop |
| 5. | Vial Jam | Machine Stop |
| 6. | No rubber stopper | Machine Stop |
| 7. | Low Air | Machine Stop |
| 8. | Vial Jam at Outfeed | Machine Stop |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | CO |
|--|---|----------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 14 of 16 | |

9 UTILITY SPECIFICATIONS

• The electrical specifications where the machine is stored should be within the ranges specified by the hardware manufacturer.

| PARAMETER | DESCRIPTION |
|-------------------------|---|
| ELECTRICAL POWER SUPPLY | 380 V, 3 PHASES, 50HZ 4 CORE, 4 SQ MM AND EARTHING |
| POWER REQUIREMENT | 4 κW |
| AIR REQUIREMENT | 4 то 6 KG/См2 ат 11 CFM |

| | Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | ASTAGENE BIOPHARMACEUTICALS CO | | |
|----------------|--|---|--------------------------------|---|--|
| | Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 | |
| Document Title | | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 15 of 16 | | |

10 REPORT APPROVAL

This Design Qualification Protocol has been completed, reviewed, and found acceptable. Any deviations found have been reconciled and this system operates and performs properly in accordance with design specifications. This Design Qualification Protocol has been reviewed and approved by the individuals listed below.

| Approval Table (Client) | | | | | |
|-------------------------|--|--|--|-----------|--|
| SIGNS FOR ROLE NAME | | | | SIGNATURE | |
| REVIEWED BY | | | | | |
| APPROVED BY | | | | | |

| Conclusion For Design Qualification: |
|--------------------------------------|
| |
| |
| |
| |
| |
| |

| Fabtech Technologies International Ltd Life Engineering | M/s. RASTAGENE BIOPHARMAG | CEUTICALS | СО |
|--|---|----------------|----|
| Document | FTIL-RASTAGENE-20-DQ-VD-002 | Revision | 0 |
| Document Title | DESIGN QUALIFICATION FOR AUTOMATIC TWO HEAD INJECTABLE LIQUID FILLING MACHINE | Page: 16 of 16 | |

11 ABBREVIATION

The following abbreviations are used in this document.

| ABBREVIATION | DEFINITIONS |
|--------------|-------------------------------------|
| cGMP | Current Good Manufacturing Practice |
| GEP | Good Engineering Practices |
| VFD | Variable Frequency Drive |
| Hz | Hertz |
| KW | Kilo Watt |
| НР | Horsepower |
| SS | Stainless Steel |
| RPM | Revolution per minute |
| MOC | Material of Construction |
| No | Number |
| IQ | Installation Qualification Protocol |
| OQ | Operation Qualification Protocol |
| FAT | Factory Acceptance Test |
| MM | Millimeter |
| ML | Milliliter |
| PW | Purified Water |
| WFI | Water for Injection |