


How does a magnetically driven agitator work?



:II. General Specifications

L&B

Magnetic coupled stirrer for pharmaceutical use. Impeller of aseptic design aimed for CIP / SIP. Journal bearing with hydrodynamic design for . pharmaceutical process
Impeller, open design. Does not require to be submerged during cleaning
.Ra<0.38µm



Material:

Impeller and welding plate AISI 316 L
Bearing: tungsten carbide
Male bearing O-rings: FPM (e.g. Viton®)
Male bearing Gasket: PTFE
Polish: Impeller and welding plate Ra<0.38µm, Optional electro polishing.

- Standard documentation in English, one set.
-